



# LILY Under Water Chest Drainage System

✓ Safety ✓ Simple ✓ Friendly

GMP, ISO 13485, ISO 14001, ISO 45001, ISO 17025 Certifications



LILY MEDICAL

## Indications

The underwater chest drainage bottle is designed to create the water seal which permits drainage of air and fluid from the pleural space but does not allow air to reenter the chest. The bottle can be used as single collection bottle with water seal or can be use at two-bottle or three-bottle system as water seal bottle.

## Features



- The bottle is transparent and durable, with clear capacity markings for easy recording and observation.
- The liquid container features accurate and clear graduation marks, with values that are easy to read.
  - 2.1. A horizontal line at the 300cc mark indicates the "UNDERWATER SEAL LINE".
  - 2.2. A line at the 1,500cc mark indicates "DO NOT FILL ABOVE THIS LINE".
- The water seal tube is located 2.5–3 cm below the 300cc marked line.
- Equipped with a check valve to effectively prevent backflow of waste liquid.
- Available in 1800ml capacities / with check valve / filter / water seal tube / white cap / sterile packaging.

## Performance

Device	Product Series	Model/Identifier	Specification
LILY Under Water Chest Drainage System	Semi-Rigid Connection Tube	CSB-013	Tubing: 7.8x12.7x1500 mm
			Bed Sheet Clip: 25.4x50.11 mm
			5 in 1 Universal Connector: ID 14 mm
			Rubber Band: ID 60 mm
	Under Water Chest Drainage Bottle full set	CSB-180	Chest Bottle: 1800 ml
			Chest Bottle: 1800 ml
			Semi Rigid Connection Tube: 7.8x12.7x1800 mm
			Spiral Male/Male Connector: ID: 7 mm
			5 in 1 Universal Connector: ID: 8 mm
			Bed Sheet Clip: 25.2x48.1 mm



GMP Certificate



ISO 13485:2016



LILY MEDICAL

TEL +886-37-612-625  
 FAX +886-37-612-720  
 E-MAIL [customer.service@lily-medical.com](mailto:customer.service@lily-medical.com)  
 WEB [www.lily-medical.com](http://www.lily-medical.com)  
 ADD NO.28-2, Shun Jeau Diann, Chunan Town.  
 Miaoli County, 35056. Taiwan R.O.C

