



LILY Suction Bag and Accessories

- ✓ Safety
- ✓ Simple
- ✓ Friendly

GMP, ISO 13485, ISO 14001, ISO 45001, ISO 17025 Certifications



LILY MEDICAL

Features



- 1 Serial Tube
- 2 Pour Spout
- 3 Patient Port
- 4 Handle
- 5 Vacuum Port



- 1 Unidirectional Valve



- 1 Pour Spout
- 2 Patient Port
- 3 Handle



- 1 Hydrophobic Suction Filter
 - >> Overflow protection
 - >> High air flow rate
 - >> Super Bacterial Filtration Efficiency (ASTM F2101-07) BFE 99.9%
- 2 Unidirectional Valve
- 3 VACUUM Port
Inline with negative pressure system to provide suction in the bag

Specifications

Model No.	Description
Suction Bag	
CSB-G10C5	1.0L, w/ Hydrophobic filter, 3 Holes
CSB-G15C5	1.5L, w/ Hydrophobic filter, 3 Holes
CSB-G20C5	2.0L, w/ Hydrophobic filter, 3 Holes
CSB-W20C5	2.0L, w/o Hydrophobic filter, 3 Holes
CSB-G105	1.0L, w/ Hydrophobic filter, 2 Holes
CSB-G155	1.5L, w/ Hydrophobic filter, 2 Holes
CSB-G205	2.0L, w/ Hydrophobic filter, 2 Holes
CSB-W205	2.0L, w/o Hydrophobic filter, 2 Holes
Canister	
CSB-C101	1.0L, Green
CSB-C151	1.5L, Green
CSB-C201	2.0L, Green or White
Trolley	
CSB-C20H12	2L, 1 gauge, 2-canister, H:75cm
CSB-C20H14	2L, 1 gauge, 4-canister, H:75cm
CSB-C20H24	2L, 2 gauge, 4-canister, H:75cm
CSB-C20L12	2L, 1 gauge, 2-canister, H:55cm
CSB-C20L14	2L, 1 gauge, 4-canister, H:55cm
CSB-C20L24	2L, 2 gauge, 4-canister, H:55cm
Connection Tube	
CSC-XXX	45cm, 180cm, 200cm, 230cm, 360cm Customized length



Suction Bag




Canister



Trolley

Accessories



CSB-A01
Regulator Module



CSB-A02
Suction Gauge,
0~700mmHg



CSB-A03
Shut Off Valve




CSB-A04
Wall Mount,
Green or White




Connection
Tube

Suction Liner Linkage

Single liner



Serial Connection



Flexible

Bag and canister can be used individually or in a serial connection up to 4 bottles

Safe

Completely sealed for endurance

User-friendly

Error-free, using color (white or green) for easy identification

Canister

Graduation lines are easy to read and minimum graduation is 50 ml



GMP Certificate



ISO 13485:2016



LILY MEDICAL

TEL +886-37-612-625

FAX +886-37-612-720

E-MAIL customer.service@lily-medical.com

WEB www.lily-medical.com

ADD NO.28-2, Shun Jau Diann, Chunan Town.

Miaoli County, 35056. Taiwan R.O.C

