

Airway Management

Laryngeal Mask Airway



Rapid Insertion, Single Patient Use and Sterile

The curved airway tube design for easy insertion and sterilized package for ease of use. A comprehensive range of single use laryngeal mask airways suitable for anaesthesia and emergency medicine.




LILY MEDICAL

Indications

- The supraglottic device (SAD) is designed blindly intubation (with a success rate of between 60% and 99%) and inserted into the hypopharynx to seal the opening of the glottis for ventilation.
- Help patients to oxygenate and ventilate in emergency or highly contagious situations, to reduce the risk of infection of medical staff.
- Widely used in various clinical indications, its versatility and ease of use are especially valuable for nursing staff engaged in anesthesia, resuscitation and intensive care.
- Can be used as a catheter to quickly and non-invasively enter the trachea of premature infants, and can avoid respiratory and extra-respiratory complications (pressure and intraventricular hemorrhage) associated with tracheal intubation.

Product Features

- Pre-inflated and soft cuff for improved fit and effective seal
- Fixation hooks help maintain proper cuff depth
- Groove design to avoid epiglottis blocking the airway
- Curve LMA with soft bite block to protect patient's teeth
- Curved design for easy insertion



Model Number

code #	Description	Max. Cuff Volume	size	Pt Weight	Endotracheal tube
CAR-S10S	Standard # 1/Neonate	< 4 ml	1	under 5 kg	
CAR-S15S	Standard # 1.5/Infant	< 7 ml	1.5	5-10 kg	
CAR-S20S	Standard # 2/Small pediatric	< 10 ml	2	10-20 kg	
CAR-S25S	Standard # 2.5/Large pediatric	< 14 ml	2.5	20-30 kg	
CAR-L30SG	90 degree curve # 3/Small adult/with Bite Block	< 20 ml	3	30-50 kg	7mm
CAR-L40SG	90 degree curve # 4/Medium adult/with Bite Block	< 30 ml	4	50-70 kg	7.5mm
CAR-L50SG	90 degree curve # 5/Large adult/with Bite Block	< 40 ml	5	70-100 kg	8mm



LILY MEDICAL CORPORATION

TEL +886-37-612-625
 FAX +886-37-612-720
 E-MAIL customer.service@lily-medical.com
 WEB www.lily-medical.com
 ADD NO. 28-2, Shun Jeau Diann. Chunan Town.
 Miaoli County, 35056 . Taiwan R.O.C.

