



Needle Free

Injection Site




3-way Stopcocks



Extension Tube



Infusion Set



Volume Metric Administration Set



Lily Medical Corporation

• 28-2, Shun Jiau Dams, Chungan Town, Ningtai 35056 Taiwan
• TEL: 886 37 612625
• FAX: 886 37 612720
• Web: www.lily-medical.com
• E-mail: customer.service@lily-medical.com





Needle Free



Needle Free

Lily Medical offers the needle free devices, the safe way to access I.V. lines.
To reduce the danger and risk of needlestick injurers.

Features

- Easy to use.
- Easy to clean.
- Easy to inspect visually.

● Injection Site

● 3-way Stopcocks

● Extension Tube

● Infusion Set

● Volume Metric Administration Set

● Customized Infusion Set with Lily Needle Free Sites



Lily Medical Corporation

- 28-2, Shun Jiau Dians, Chunan Town, Miaoli 35056 Taiwan
- TEL: 886 37 612625
- FAX: 886 37 612720
- Web: www.lily-medical.com
- E-mail: customer.service@lily-medical.com





Infusion

Volume Metric Administration Set



Y injection site



Air vented spike



Microdrip adapter



Needle free Y connector



Assembled needle free burette upper lid

Spec

- Available burette volume: 120/ 150ml
- Spike: \approx 20 drops/ml
- Microdrip adapter: \approx 60 drops/ml
- Tubing length: 180 \pm 5cm

Features

- Spike is sharp and easy to puncture.
- Connector fixes to each type of spike and connector.
- Roller clamp is superior control of the flow rate.
- With safety flow valve.
- The end of set can be optional to luer lock or slip.
- Standard sterilized package.

P/N	Spec
1102001	120ml/20 drops/ Needle free Y connector/ Male luer lock/Microdrip adapter/ 160cm
1102018	120ml/20 drops/ Y connector/ Male luer slip/ Microdrip adapter/ 160cm
1102023	120ml/20 drops/ Microdrip adapter/ Needle free Y connector/ Male luer slip/ 160cm
1102025	120ml/20 drops/ Y injection site/ Microdrip adapter/ Flash bulb connector/ 160cm
1102026	120ml/20 drops/ Microdrip adapter/ Needle free Y connector/ Male luer lock/ 160cm
1102029	150ml/20 drops/ Microdrip adapter/ Y connector/ Male luer lock/ 165cm

Lily Medical Corporation



- 28-2, Shun Jiau Diann, Chunan Town, Miaoli 35056 Taiwan
- TEL: 886 37 612625
- FAX: 886 37 612730
- Web: www.lily-medical.com
- E-mail: customer.service@lily-medical.com





Infusion

Infusion Set



Y injection site



Microdrip adapter



Male luer slip



Needle free Y connector



Air vented spike

Spec

- Spike: ≈ 20 drops/ml
- Chamber: ≈ 20 drops/ml
- Total length: ≈ 180cm

Features

- Spike is sharp and easy to puncture.
- Connector fixes to each type of spike and connector.
- Roller clamp is superior control of the flow rate.
- The end of set can be optional to luer lock or slip.
- Standard sterilized package.

P/N.	Spec
1103033	3.1x4.1x1650mm/ Needle free Y connector/ Male luer slip
1103066	3.1x4.1x1650mm/ Needle free Y connector/ Microdrip adapter/Male luer slip
1103010	3.1x4.1x1390mm/ 2xNeedle free Y connector/ Male luer lock(C-spike)
1103011	3.1x4.1x1500mm/ Needle free Y connector/ Male luer lock(C-spike)
1103012	3.1x4.1x1500mm/ Needle free Y connector/ Microdrip adapter/ Male luer lock(C-spike)
1103013	3.1x4.1x1390mm/ 2xNeedle free Y connector/ Microdrip adapter/ Male luer lock(C-spike)

Lily Medical Corporation

- 28-2, Shun Jiau Diann, Chunan Town, Miaoli 35056 Taiwan
- TEL: 886 37 612625
- FAX: 886 37 612720
- Web: www.lily-medical.com
- E-mail: customer.service@lily-medical.com





Infusion

Extension Tube



P/N.	Spec
1106004	Two ways/ 1.67mmx3.5mmx80mm/ Male luer slip/ Female luer slip
1106006	2.5mmx4.1mmx110mm/ Needle free Y connector/ Male luer lock
1106008	1.2mmx2.4mmx110mm/ Needle free Y connector/ Male luer lock
1106009	3.1mmx4.1mmx190mm/ 0.22 μ Filter/ Needle free Y connector/ Male luer lock
1106011	Four ways/ 1.67mmx3.5mmx80mm/ Male luer lock/ Female luer slip
1106013	3.1mmx4.1mmx250mm/1.2 μ Filter/ Needle free Y connector/ Male luer lock
1106014	1.2mmx2.4mmx60mm/ 2xNeedle free injection site/ 2-port connector/ Male luer lock
1106015	3.1mmx4.1mmx75mm/ 2xNeedle free injection site/ 2-port connector/ Male luer lock
1106018	3.1mmx4.2mmx30mm/ Vented Spike/ Check Valve/ Needle free Injection site
1106059	1.27mmx2.4mmx2.4mmx100mm/ Male luer lock/ Female luer slip
1106056	1.27mmx2.4mmx170mm/ Needle Free injection site/ Male luer lock
1106060	1.27mmx2.4mmx170mm/ Needle Free injection site/ Male luer slip

Lily Medical Corporation

- 28-2, Shun Jiau Diann, Chunan Town, Miaoli 35056 Taiwan
- TEL: 886 37 612625
- FAX: 886 37 612720
- Web: www.lily-medical.com
- E-mail: customer.service@lily-medical.com





Infusion

T-connector



Features

- Luer lock design; usable for a vein injection patient.
- Tube won't be out of shape after bending, and the connector can be fixed tightly with each other.

P/N	Spec
1106050	1.27x2.4x150mm/ T connector/ Female luer slip/ slide clamp
1106051	1.27x2.4x320mm/ T connector/ Female luer slip/ slide clamp
1106061	1.27x2.4x170mm/ T connector/ Needle free injection site/ Female luer slip/ slide clamp

Lily Medical Corporation

- 28-2, Shun Jiau Diann, Chunan Town, Miaoli 35056 Taiwan
- TEL: 886 37 612625
- FAX: 886 37 612720
- Web: www.lily-medical.com
- E-mail: customer.service@lily-medical.com





Lily Medical

Infusion

CVP Manometer



Optional for
3-Way stopcock



3-way stopcocks
180° turning
white handle



3-way stopcocks
360° turning
blue handle

Features

- -3 to -3.8cm H₂O.
- With a ball in the tubing, easy to view.
- PCT material tubing is with characters of hardness, transparence, easy to do the after processing, eco-friendly, no odor, non-toxic, and will not cause pollution when burning.
- With a stand.
- 180° or 360° optional 3-way stopcocks.

P/N.	Spec
1107001	-3 to 35cm/ 125cm tube/3-way stopcocks with blue handle/Y injection site/ Male luer slip/ Female luer slip/Flying ring/Red holder/White stand
1107002	-3 to 35cm/ 125cm tube/3-way stopcocks with white handle/Y injection site/ Male luer lock /Female luer slip /Flying ring/Red holder/White stand

Lily Medical Corporation

- 28-2, Shun Jiau Diann, Chunian Town, Miaoli 35056 Taiwan
- TEL: 886 37 612625
- FAX: 886 37 612720
- Web: www.lily-medical.com
- E-mail: customer.service@lily-medical.com





Lily Medical

Suction

Collection bottle



- P/N: 1104007 / 1204001
- Single use collection bottle

Purpose

- Design for collection body fluids from the patients.

Feature

1. Clear bottle to easy inspect visually.
2. Moveable lid.
3. Available sizes: 1750 or 3000ml capacity
4. Optional filter or back flow valve.



- P/N: 1104001 / 1104002
- Chest drainage bottle

Purpose

- Normal used after cardiac and thoracic surgery.

Feature

1. Clear bottle to easy inspect visually.
2. Water-seal one-way design permits air and fluid to be removed while preventing backflow into the chest.
3. Filter/back flow valve/under water sealed tube/white lid/1800 or 3000 ml



- P/N: 1104003
- Woung drainage bottle

Purpose

- Design for collection body fluids from the patients; coordinate to use with chest drainage bottle.

Feature

1. Clear bottle to easy inspect visually.
2. Filter protects filtration and avoid overflow.
3. 1800ml/filter/backflow valve/blue lid



- P/N: 1104004
- Reusable collection bottle

Purpose

- The bottle is used to collect waste liquid or mucus from medical suction.

Features

1. Construer can be withstand at 135°C or 2% Lysol.
2. The bottle is limp and solid; scale marks clear, it is easy to record and observe.
3. The sealed bottle lids can drainage immediate to save time and be safe.
4. Non-return valve effectively prevents reflux retrograde infection.
5. 1500ml/Backflow valve/Filter/Red lid

Lily Medical Corporation

- 28-2, Shun Jiau Diann, Chunan Town, Miaoli 35056 Taiwan
- TEL: 886 37 612625
- FAX: 886 37 612720
- Web: www.lily-medical.com
- E-mail: customer.service@lily-medical.com





Lily Medical

Suction

Collection tube & Yankauer suction set



Features (tubing)

- Made of PVC, not easy to be break.
- Tube is transparent, easy to view.
- Not easy to be bended by waving design.

P/N.	Spec
1110002	7.0x10.4x450mm
1110011	7.0x10.4x1800mm
1110012	7.0x10.4x2300mm
1110013	7.0x10.4x3600mm
1210010	7.0x10.4x1800mm
1210005	7.0x10.4x2300mm
1210007	7.0x10.4x3600mm

Lily Medical Corporation

- 28-2, Shun Jiau Diann, Chunan Town, Miaoli 35056, Taiwan
- TEL: 886 37 612625
- FAX: 886 37 612720
- Web: www.lily-medical.com
- E-mail: customer_service@lily-medical.com





Lily Medical

Suction

Suction Catheter



Purpose

- Remove the secretions from the trachea and bronchial.
- Keep the airway clear and prevent plugging.

Features

- Clear tubing enables to view the secretions easily.
- Standard color-coded controller for instant identification.
- Including one glove.
- The top of tube is spiral type which can be linked to kinds of tube machine.
- Thumb controller to fit the best suctioning needs.

P/N.	Spec
1201034	06Fr/1.19x2.0x470mm/with glove
1201006	08Fr/1.78x2.7x470mm/with glove
1201008	10Fr/2.19x3.3x470mm/with glove
1201010	12Fr/2.68x4.0x500mm/with glove
1201013	14Fr/3.18x4.7x470mm/with glove
1201016	16Fr/3.72x5.0x470mm/with glove

Lily Medical Corporation

- 28-2, Shun Jiau Diann, Chunan Town, Miaoli 35056 Taiwan
- TEL: 886 37 612625
- FAX: 886 37 612720
- Web: www.lily-medical.com
- E-mail: customer.service@lily-medical.com





Lily Medical

Suction

Mucus Extractor



Purpose

- For cultivating the bacteria of sputum from the patient who operate the tracheostomy or endotracheal intubation.
- For examining properties or abnormal cell in the sputum.

Features

- With aspiration controller.
- Clear scale on the limp bottle for easy observation.
- Smooth surface minimize the outer friction.

P/N.	Spec
1101004	06Fr/1.19x2.0x470mm
1101007	08Fr/1.78x2.7x470mm
1101011	10Fr/2.19x3.3x470mm
1101015	12Fr/2.68x4.0x500mm
1101019	14Fr/3.18x4.7x470mm
1101023	16Fr/3.72x5.0x470mm

Lily Medical Corporation

- 28-2, Shun Jiau Diann, Chunan Town, Miaoli 35056, Taiwan
- TEL: 886 37 612625
- FAX: 886 37 612720
- Web: www.lily-medical.com
- E-mail: customer.service@lily-medical.com





Lily Medical

Drainage Bag

Urine Bag



T-type urine bag
P/N: 1109007

Features

- T type slucing valve of bottom outlet.
- Conical inlet connector with cap.
- Soft and kink resistant tubing.
- Single package with sterilized.

Spec: 2000ml/100cm tubing/T-type slucing valve



Slucing valve



Standard urine bag
P/N: 1109008

Features

- Smooth and kink resistant tubing.
- Easy handled by single hand of bottom outlet.
- Clear measuring volume
- Air vent secures constantly airflow through system.
- Optional sampling port.
- Single package with sterilized.

Spec: 2000ml/120cm tubing/sampling port



Sampling Port



Air vent

Lily Medical Corporation

- 2B-2, Shun Jiau Diann, Chunan Town, Miaoli 35056 Taiwan
- TEL: 886 37 612625
- FAX: 886 37 612720
- Web: www.lily-medical.com
- E-mail: customer.service@lily-medical.com

